
PULSE OXIMETER

Baby Oxygen Monitor

User Manual (EN)

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Manuale d'uso (IT)

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Contents

User Manual.....	English 1-24
Benutzerhandbuch.....	Deutsch 25-50
Manuale d'uso.....	Italiano 50-75
Manual de usuario.....	Español 75-99
Manuel de l'utilisateur.....	Français 99-124

Contents

1 Tips for Safe Sleep	2
2 Warnings and Cautions	2
3 Name and Model	4
4 Unpacking	4
5 Overview	5
6 Intended Use and Contraindications	6
7 Using the device	6
8 Maintenance	12
9 Troubleshooting	12
10 Specifications	14
11 Symbols	15
12 FCC Statement	17
13 Appendix: Electromagnetic compatibility	18

1 Tips for Safe Sleep

- Room-share instead of co-sleeping.
- Place your baby alone and on their back to sleep.
- Use a firm sleep surface. Keep soft objects out of crib.
- Use a crib that meets safety standards. Do not use cradle for sleep.
- Don't use the device as an excuse for unsafe sleep.

2 Warnings and Cautions

- DO NOT twist the sensor and the wrap or apply excessive force on them.
- Do not use this device during MRI examination.
- Do not store the device in the following locations: locations in which the device is exposed to direct sunlight, high temperatures or levels of moisture, or heavy contamination; locations near to sources of water or fire; or locations that are subject to strong electromagnetic influences.
- Do not use the device in a combustible environment.
- Never submerge the device in water or other liquids.
- Do not clean the device with acetone or other volatile solutions.
- Do not drop this device or subject it to strong impact.
- The device and accessories are provided non-sterile.
- Do not place this device in pressure vessels or gas sterilization device.
- Do not dismantle the device, as this could cause damage or malfunctions or impede the operation of the device.
- Consult your doctor immediately if you experience symptoms that could indicate acute disease.

- Do not self-diagnose or self-medicate on the basis of this device without consulting your doctor. In particular, do not start taking any new medication or change the type and/or dosage of any existing medication without prior approval.
- Use only cables, sensors and other accessories specified in this manual.
- Prolonged continuous monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering or burns.
- Do not disassemble or modify the device without authorization of the manufacturer, otherwise it may cause machine malfunction or affect the normal operation of the device.
- The device shall only be maintained by qualified professionals.
- The biocompatibility testing has been performed on the materials in contact with the person in accordance with ISO10993.
- Do not place the SpO₂ probe on a finger or toe with edema or fragile tissue.
- Check the SpO₂ sensor and cable before use. Do not use a damaged SpO₂ sensor.
- Check the SpO₂ sensor application site every 6-8 hours to determine the positioning of the sensor and the circulation and skin sensitivity of the patient. Patient sensitivity varies depending on medical status or skin condition. For patients with poor peripheral blood circulation or sensitive skin, inspect the sensor site more frequently.
- The functional tester cannot be used to assess the accuracy of the SpO₂ sensor or a device.
- The device has no alarm system.

- The local laws and regulations should be followed when disposing of the device and accessories.
- The battery inside this product cannot be replaced
- The manufacturer shall provide the service personnel with circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist service personnel to repair the device.
- The patient is the expected user.

3 Name and Model

Product name: Pulse Oximeter


Model: PO5, PO5-1, PO5-2, PO5-3

Model	Enclosure color	Intended users	Measuring part
PO5	Blue	Adults, pediatric and infants	Finger/Foot
PO5-1	Pink		
PO5-2	White		
PO5-3	Green		

4 Unpacking

- Device (main unit)
- Sensor cable
- Soft wraps
- User Manual
- Quick guide
- Data/Charging Cable

5 Overview

1. Device

Displays device status and measurement results.

2. SpO2 sensor

3. Soft wraps

Secure the sensor and main unit to the measuring site.

4. USB data/charging cable

Export measurement data or charge the device.

5. Side button

- Press and hold 2 seconds to power on/off.
- Short press to wake up the screen or switch display.

6. Display screen

6 Intended Use and Contraindications

6.1 Intended Use

This Pulse Oximeter is intended to be used for measuring, displaying and storing of the pulse oxygen saturation (SpO₂), pulse rate of adults, pediatric and infants in home or healthcare facilities environment.

6.2 Contraindications

No contraindications.

7 Using the device

7.1 Charging

Charge the battery before use.

Connect the device to a computer or USB charging adapter with the supplied cable.

The device will power off automatically when fully charged.

Note: The device cannot be used during charging, and if choosing a third party charging adaptor (Class II), select one that complies with IEC62368 or IEC60601-1.

7.2 Power on/off

Power on:

Press the side button for 2 seconds.


Power off:

Three ways:

- Press the side button for 2 seconds.
- Unplug the sensor cable, it will power off automatically.
- With the sensor cable connected, it will power off automatically if not in use for 10 minutes.

7.3 Placement Guide before working

The device and sensor cable are packaged with four split wraps (three foot wraps of different sizes and a leg wrap).

To ensure that the measurement is performed normally, please follow the steps below to wear the device.


1) Select a foot wrap

Choose the proper foot wrap based on the baby's weight.

Size	Weight
S	<3.9kg (8.6 lbs)
M	3.9-8.5kg (8.6~18.7 lbs)
L	8.5-13.9kg (18.7~30.6 lbs)


2) Place the sensor

Plug the sensor into the wrap as shown below.


3) Secure the sensor

- A. Place the foot wrap on the right foot with the sensor part covering the outer side of the foot. The sensor cable should be oriented to the baby's leg.


- B. The light emitter and the light detector should be opposite to each other on the two sides of the foot. Attach the two Velcro on the highly elastic wrap securely with no gap between the sensor and the baby's skin.

Note: The light emitter should be on the top while the light detector at the bottom.

4) Attach the monitor

Plug the sensor cable into the device. Then attach the device in a proper position on the right leg to ensure the cable is not too tight or too loose.


7.4 Start working


After wearing the device correctly, please follow the steps below to measure oxygen level and pulse rate.

- 1) Press the side button for 2 seconds to turn on the device. The screen displays time and battery. There should be red light inside the soft wrap.
- 2) Press the side button to set the oxygen threshold.

After a few seconds, readings will be displayed on the screen.

If the screen displays “poor signal”, please check if the sensor has been positioned and secured properly.

When the above symbol displays on device screen, it indicates the readings are unavailable right now.

When at rest, there will be readings in a few seconds. If you cannot get readings for a long time, it may be caused by:

- Improper sensor placement
- Frequent movement
- Foot is too cold

Notice:

- If the working time is less than 2 minutes, the data will not be saved.
- Motions may cause readings unavailable, which will recover very soon when the baby is at rest.
- Please avoid strong ambient light condition.


SpO2 measurement principle:

The Pulse Oximeter is a lightweight, portable health oximeter for use in the home or in healthcare facilities. SpO2 measurement technology is based on developed photoelectron method, the circuit design and calculation software were developed by Shenzhen Viatom Technology Co., Ltd. According to Lambert-Beer law, when the hand or foot is placed under a light source, the difference between the transmitted light intensity and the incident light intensity measured is the absorbed light energy, which contains

information about the proportion of oxygenated hemoglobin in the user's blood. The MCU receives the pulse signal, gets the frequency signal by counting, processes its digital signal, and finally gets the measured SpO2 value. The PR is averagely calculated by above peak intervals of PR waveform.

7.5 Stop working & sync data

Unplug the sensor cable, the countdown on the screen will begin.

Stop ? 10

During the countdown, if you plug the sensor cable again, the record will be resumed.

After the countdown, the data will be ready for sync.

Notice:

- If the working time is less than 2 minutes, there will be no data saved
- If the device power off during working, the data should have been saved.

Sync data:

- After the countdown, open the Bluetooth or through the USB port to sync data.
- You can also sync data later. The built-in memory can store 4 sessions. The oldest will be overwritten by the 5th.

Notice: The built-in memory can store 4 sessions, up to 10 hours for each. It will start the next session automatically after 10 hours monitoring.

7.6 Screen wake up / change display

During measurement, the screen will go off automatically for saving power; you can press the side button to wake up the screen or change display.

7.7 Bluetooth Connection

The device Bluetooth will be enabled automatically after it's turned on.

To establish a Bluetooth connection,

- 1) Power on the device.
- 2) Make sure the phone Bluetooth is enabled.

8 Maintenance

8.1 Time & Date

After Bluetooth connection, the device time will be synced automatically.

8.2 Cleaning

Clean the wrap with water. Water temperature should be lower than 104°F (40°C)

Use a soft cloth moistened with water or alcohol to clean the device surface.

Note: The device is a non-sterile medical device and does not contain any sterile or degradable component thus the device is not subject to the shelf-life requirements.

8.3 Recycle

Relevant wastes, residues, and end-of-life equipment and accessories should not be discarded arbitrarily, and should comply with local regulations. When you intend to dispose of this device, you must send it to an appropriate facility for recovery and recycling.

9 Troubleshooting

Problem	Possible Cause	Possible Solution
No readings	the sensor is not positioned and	check if the sensor has been positioned and

on screen	secured properly	secured properly
	Baby is moving frequently	It will recover when at rest
Device does not turn on or no response	Battery may be low.	Charge battery and try again.
	Device might be damaged.	Please contact the seller or the after-service email.
The device screen displays “Error 1”.	Errors occur during data analysis.	Connect the power supply and plug in the charging cable, press and hold the touch key for 3s to reset the hardware.
The device screen displays “Error 4”.	Sensor error.	Please re-plug the sensor.
Bluetooth connection error	The Bluetooth of your phone is off.	Turn on the Bluetooth in the phone.
	The device Bluetooth is off.	Turn on the device
	For Android, Bluetooth cannot work without location permission	Allow location access


10 Specifications

Classification		
Regulation or Directives	EU 2017/745	
	RED Directive 2014/53/EU	
Environmental	Operating	Storage
Temperature	5 to 40°C	-25 to 55°C
Relative humidity (noncondensing)	10% to 80%	10% to 93%
Barometric	700 to 1060hPa	700 to 1060hPa
Protection against electric shock	Internally powered equipment	
Degree protection against electrical shock	Type BF	
Electro-magnetic compatibility	Group I, Class B	
Degree of dust & water resistance	IP22	
Weight	31g (main unit with band and sensor)	
Size	49mm×27mm×14 mm (main unit)	
Battery	3.7 Vd.c., Rechargeable Li-ion Polymer Battery	
Charge time	2-3 hours	
Battery life	16 hours for typical use	
Wireless	Bluetooth 4.2 BLE	

Oxygen level range	0% to 100%
Oxygen level Accuracy (Arms)	70%-100%: $\pm 2\%$, 0%-69%: Undefined
Pulse rate range	20 to 250 bpm
Pulse rate accuracy	± 2 bpm or $\pm 2\%$, whichever is greater
Recorded parameters	Oxygen level, heart rate
Wavelength/Max emission power	660nm/940nm, 0.8mW/1.2mW
Data storage	4 sessions, up to 10 hours for each
Expected service life	5 years
Bluetooth	
Frequency Range	2.402 – 2.480 GHz
Modulation type	GFSK modulation
Max RF power	8 dB

11 Symbols

Symbol	Description
	Manufacturer
	Date of manufacture
	Distributor
	Medical device

	Serial number
	Unique device identifier
	Indicates a medical device that is not to be disposed of as unsorted municipal waste.
	Follow Instructions for Use.
	Type BF Applied Part
	No alarm system
	MRI unsafe. Presents hazards in all MR environments as device contains strongly ferromagnetic materials.
IP22	Resistant to liquid ingress
CE 0197	CE marking
	Authorized representative in the European community
UK CA	UKCA marking
	Authorized Representative in the United Kingdom
FC	This product complies with the rules and regulations of the Federal Communication Commission.
	Non-ionizing radiation
	This product complies with verpackG.

	<p>Our products and packaging can be recycled, don't throw them away!</p> <p>Find where to drop them off on the www.quefairedemesdechets.fr site (Only applicable for French market).</p>
	<p>Temperature limitation</p>
	<p>Humidity limitation</p>
	<p>Atmospheric pressure limitation</p>

12 FCC Statement

FCC Warning:

FCC ID: 2AD XK-4728

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and

used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

13 Appendix: Electromagnetic compatibility

The device meets the requirements of IEC 60601-1-2.

Warnings and Cautions

This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.

The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

Guidance and manufacturer's declaration

– electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	n.a.	
Voltage fluctuations/ flicker emissions IEC 61000-3-3		

Guidance and manufacturer’s declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC61000-4-6	3Vrms 150kHz to 80MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC61000-4-3	10V/m 80MHz to 2.7GHz	10V/m	

		<p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p>80MHz to 800MHz</p> $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>800MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 - At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2-These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>		
<p>a The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4</p>		

MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/	n.a.	n.a.

	output lines		
Surge IEC61000-4-5	± 1 kV line to line ±2 kV line to earth	n.a.	n.a.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT 0,5cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°, 0% UT 1cycle and 70% UT 25/30 cycles Single phase: at 0°	n.a.	n.a.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m, 50/60Hz	30A/m, 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the AC mains voltage prior to application of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and the A&D unit

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of	Separation distance according to frequency of transmitter (m)
-------------------------------	---

transmitter (W)	150kHz to 80MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80MHz to 800MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800MHz to 2.7GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.04	0.07
0.1	0.37	0.12	0.23
1	1.17	0.35	0.70
10	3.70	1.11	2.22
100	11.70	3.50	7.00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recommended separation distances between RF wireless communications equipment

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Frequency MHz	Maximum Power W	Distance	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
385	1.8	0.3	27	27	RF wireless

450	2	0.3	28	28	communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
710	0.2	0.3	9	9	
745					
780	2	0.3	28	28	Recommended separation distance
810					
870					
930	2	0.3	28	28	$E = \frac{6}{d} \sqrt{P}$
1720					
1845					
1970	2	0.3	28	28	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
2450					
5240					
5500	0.2	0.3	9	9	
5785					
5785					

NOTE 1 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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MedNet EC-REP GmbH
MediMap Ltd
Australia Sponsor: SHARE INFO PTY LTD